

# NEURO ATTACHMENT



**Three Pin Head Neuro Attachment**



**Horse Shoe Shaped Neuro Attachment**



## Objectives

We offer a supreme range of Neuro Attachment, which is suitable for brain surgery. This Neuro attachment has been developed & manufactured with innovative quality to meet the demanding needs of Neurosurgeons.

## Key Specifications

### Three Pin Head Neuro Attachment: -

- The skull clamp has the additional feature of a swiveling C-ARM.
- This swiveling capacity greatly facilitates positioning the pins on the head around nerves, vessels, and thin bone in the temporal region.
- For better result, the arm is locked by control knob and the width is adjusted by a ratchet mechanism.
- The three cone-shaped pins available in adult and child sizes are interchangeable and may be removed for sterilization.

### Horse Shoe Head Shaped Neuro Attachment: -

- This attachment, originally designed for transphenoidal hypophysectomy is used for all prone position procedures.
- It is attached to the Adjustable Basic Unit by means of transitional member and the Gardner Clamp Adaptor.
- It is padded with soft, pliable sponge rubber for equalized weight distribution and incorporates the new quick-change pads.
- A removable, adjustable pulley bar is included if traction is desired.



---

**For Further Information, contact:**



487/67, 2<sup>nd</sup> Floor, Glass Street,  
National Market, Peeragarhi,  
New Delhi, India – 110087

Tele: 91-11-25255415 / 416

Fax: 91-11-25255417

Mail: [info@meditech-india.com](mailto:info@meditech-india.com)

Web: [www.meditech-india.com](http://www.meditech-india.com)

**Online assistance:**



MSN: marketing\_meditech



Skype: exportmeditech