



Mini 2



Mini 1

Technical Data

	Mini 2	Mini 1/Mini M
Light Intensity @ 1m	80,000 + 80,000 Lux ±10%	80,000 + 80,000 Lux ±10%
No. of LED's	15+15	15
LED's Lifetime	>50000 Hrs.	>50000 Hrs.
Dome Diameter	400mm x 2	400mm
Light Body Type	Fire Resistant Polymer	Fire Resistant Polymer
Colour Rendering Index	93 RA	93 RA
Temp. Increase @ Head Area	2° C	2° C
Colour Temperature	4500K	4500K
Focus Adjustment	Fixed	Fixed
Brightness Control	Capacitive Touch Panel	Capacitive Touch Panel
Remote Control	Yes	Yes
Light Field Diameter	15 cm	15 cm
Colour of LED's	White & Yellow (Double Colour)	White & Yellow (Double Colour)
Power Supply	220-230V AC/50-60Hz	220-230V AC/50-60Hz



Mini M