



ELECTROMATIC O. T. TABLE MI 103





Overview

This Electromatic O.T Table is suitable for all kind of surgeries, i.e. general surgery, vascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, dental surgery, plastic surgery and others. The table is specially designed with hanging ortho attachment to support a patient in correct positions for various orthopedic procedures.

Key Structure

- The five section table top is X-ray translucent and is equipped with removable, antistatic 50 mm thick mattresses.
- The base & column of the table is made of high quality medical stainless steel (Grade - 304). Non- reflecting surface, antibacterial & is easy to clean.
- The table can provide with a top sliding of 250 mm, Dual Leg Section & Battery backup as per consumer requirement. **(Optional)**
- Hanging Orthopedic Attachment including Pelvic Support, Tibia Nailing & Femur Nailing enables the limbs to perform strain free.
- Height adjustment, Lateral Tilts, Trendelenburg, and Flex - Reflex, Chair Position and Longitudinal Sliding top are realized by Remote control.
- Head and Leg sections are detachable & interchangeable.

Technical Constraints	
Top dimension	L 2032 x W 550 mm
Height adjustment	787 mm - 1037 mm
Table Top Sliding	250 mm (optional)
Trendelenburg / Reverse	30° / 25°
Lateral tilt	20° / 20°
Kidney elevator	150 mm
Back Rest (up / down)	80° / 25°
Leg Rest (up / down)	15° / 90°
Head Rest (up / down)	20° / 60°
Power supply	24V DC
Battery Backup	2 - 3 hours

TABLE POSITION



ACCESSORIES





OPERATING TABLE FUNCTIONS

