



Lexia 2



Lexia 1

Technical Data	Lexia 2	Lexia 1/Lexia M
Light Intensity @ 1m	1,75,000 + 1,75,000 Lux ±10%	1,75,000 Lux ±10%
No. of LED's	105+105	105
LED's Lifetime	>50000 Hrs.	>50000 Hrs.
Dome Diameter	650mm x 2	650mm
Dome Type	Fire Resistant Polymer	Fire Resistant Polymer
Colour Rendering Index	93 RA	93 RA
Temp. Increase @ Surgeon Head	2° C	2° C
Colour Temperature	4000K-5000K	4000K-5000K
Focus Adjustment	Adjustable Manually	Adjustable Manually
Brightness Control	Capacitive Touch Panel	Capacitive Touch Panel
Light Field Diameter	12-30 cm	12-30 cm
Colour of LED's	White & Yellow (Double Colour)	White & Yellow (Double Colour)
Power Supply	220-230V AC/50-60Hz	220-230V AC/50-60Hz
Suspension Type	Stainless Steel	Stainless Steel

Using The Best LED Technology



Lexia M