

ORTHOPEDIC ATTACHMENT MI 112



Hip Replacement



Tibia Nailing



Femur Nailing



Objectives

Ortho Attachment allows the table to be used for the lower extremity procedures. Available at competitive prices our ortho attachment is purely made of stainless steel which resist corrosion.

Key Specifications

- This Orthopedic Extension device can be used with any operating table as it independently stands on the floor.
- Specially Design for Intra Operative X-Ray photos & image intensifier set use.
- In addition, a pelvic support plate is supported on top of the frame and is formed with an asymmetrical construction. The pelvic support plate includes a cutout area on a front corner thereof permitting additional access to the affected leg of a patient. Removable Pelvic rest & X-Ray translucent pelvic stem.
- Stirring Instruments allows Flexion/Raising of Legs & Stirring/String ball joints locked from foot end.
- Two Traction Devices on sliding spot with ball joint allowing flexions & rotations of foot with traction shoe with metallic sole.
- Padded orthopedic knee Crutches/Provides support to the legs of the patient.
- Traction Unit with supporter is so easily & collapsible that it requires a very small place.
- The plate may be flipped over to position the cutout area on the opposite side of the patient, depending on the location of the affected leg.



For Further Information, contact:



487/67, 2nd Floor, Glass Street,
National Market, Peeragarhi,
New Delhi, India – 110087

Tele: 91-11-25255415 / 416

Fax: 91-11-25255417

Mail: info@meditech-india.com

Web: www.meditechindia.in

Online assistance:



MSN: marketing_meditech |



Skype: exportmeditech