



**MEDITECH**  
INDIA PVT. LTD.

SWITCH TO   
**BETTER**  
Suction Machine




ISO  
9001:2015  
Certified  
Company

ISO  
13485:2016  
Certified  
Company



<b>Housing</b>	Powder coated body (Trolley Type)
<b>Max. Vacuum</b>	-700mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
<b>Motor</b>	<b>Oil free double piston Pump</b>
<b>Suction Capacity</b>	<b>35 Liters / Min</b>
<b>Vacuum Gauge</b>	2.5 inch, 0-760 mm Hg
<b>Jars</b>	2 x 1.5 Liters Polycarbonate Jars
<b>Tubing</b>	Non-Collapsible medical grade tubing
<b>Noise Level</b>	<30-35dB A $\pm$ 3 Almost Silent
<b>Overflow Safety</b>	Mechanical
<b>Certifications</b>	CE Certified

## VacuMatrix VM - 35






**MEDITECH**  
INDIA PVT. LTD.

over  
**30**  
YEARS

SERVICE & EXPERIENCE

Housing	Powder coated body (Trolley Type)
Max. Vacuum	-720mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
Motor	<b>Oil free double piston Pump</b>
Suction Capacity	<b>50 Liters / Min</b>
Vacuum Gauge	3 inch, 0-760 mm Hg
Jars	2 x 2 Liters Polycarbonate Jars
Tubing	Non-Collapsible medical grade tubing
Noise Level	<40-45dB A $\pm$ 3 Almost Silent
Overflow Safety	Mechanical
Certifications	CE Certified

## VacuMatrix VM - 50



A trusted Indian manufacturer and global exporter of Operating Tables, Surgical & Examination Lighting, Suction Machine, Fogger and Hospital Furniture - serving healthcare since 1996



**MEDITECH**  
INDIA PVT. LTD.

over  
**30**  
YEARS

SERVICE & EXPERIENCE

<b>Housing</b>	Powder coated body with Stainless Steel Top
<b>Max. Vacuum</b>	-710mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
<b>Suction Capacity</b>	30 -35 Liters / Min
<b>Vacuum Gauge</b>	2.5 inch , 0 -760mm Hg
<b>Jars</b>	2 x 1.5 Liters Poly Carbonate jars or Glass jars ( Optional )
<b>Tubing</b>	Non- Collapsible medical grade tubing
<b>Motor</b>	Goley motor ¼ hp
<b>Noise Level</b>	<50dB A $\pm$ 3 Almost Silent
<b>Vacuum pump</b>	Rotary vane Type
<b>Certifications</b>	CE Certified

## VacuMatrix VM - 30



A trusted Indian manufacturer and global exporter of Operating Tables, Surgical & Examination Lighting, Suction Machine, Fogger and Hospital Furniture - serving healthcare since 1996



**MEDITECH**  
INDIA PVT. LTD.

over  
**30**  
YEARS

SERVICE & EXPERIENCE

<b>Housing</b>	Powder coated body with Stainless Steel Top
<b>Max. Vacuum</b>	-730mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
<b>Suction Capacity</b>	40 -55 Liters / Min
<b>Vacuum Gauge</b>	3 inch , 0 -760mm Hg
<b>Jars</b>	2 x 2 Liters Poly Carbonate jars or Glass jars ( Optional )
<b>Tubing</b>	Non- Collapsible medical grade tubing
<b>Motor</b>	Goley motor 1/2 hp
<b>Noise Level</b>	<50dB A $\pm$ 3 Almost Silent
<b>Vacuum pump</b>	Rotary vane Type
<b>Certifications</b>	CE Certified


## VacuMatrix VM - 55



A trusted Indian manufacturer and global exporter of Operating Tables, Surgical & Examination Lighting, Suction Machine, Fogger and Hospital Furniture - serving healthcare since 1996

<b>Housing</b>	Powder coated body with Stainless Steel Top
<b>Max. Vacuum</b>	-740mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
<b>Suction Capacity</b>	70 - 75 Liters / Min
<b>Vacuum Gauge</b>	3 inch , 0 -760mm Hg
<b>Jars</b>	2 x 2 Liters Poly Carbonate jars or Glass jars ( Optional )
<b>Tubing</b>	Non- Collapsible medical grade tubing
<b>Motor</b>	Goley motor 1 hp
<b>Noise Level</b>	<50dB A $\pm$ 3 Almost Silent
<b>Vacuum pump</b>	Rotary vane Type
<b>Certifications</b>	CE Certified

## VacuMatrix VM - 75






**MEDITECH**  
INDIA PVT. LTD.

over  
**30**  
YEARS

SERVICE & EXPERIENCE

Housing	Powder coated body (Trolley Type)
Max. Vacuum	-760mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
Motor	<b>Oil free double piston Pump</b>
Suction Capacity	<b>80 Liters / Min</b>
Vacuum Gauge	3 inch, 0-760 mm Hg
Jars	2 x 2 Liters Polycarbonate Jars
Tubing	Non-Collapsible medical grade tubing
Noise Level	< 50dB A $\pm$ 3 Almost Silent
Overflow Safety	Mechanical
Certifications	CE Certified

## VacuMatrix VM - 80



VACUUM PUMP



FOOT SWITCH



COOLING FAN




PC JAR 2000 ML






**MEDITECH**  
*INDIA PVT. LTD.*

 #487/67, Glass Street, National Market,  
Peeragarhi, New Delhi - 110087, INDIA

 +91 11 40392484

 +91 9812266676

 [marketinguk@meditech-india.com](mailto:marketinguk@meditech-india.com)

 [www.mediech-india.com](http://www.mediech-india.com)

\*\* This document is intended for the exclusive use of Meditech India Pvt. Ltd. Customers, including architects or designers. Reproduction in whole or in part by any party other than a Customer is prohibited.